

Datasheet

Goat Anti-Chicken IgG (H&L) secondary antibody

Catalog Number: PAB9169

Regulation Status: For research use only (RUO)

Product Description: Goat anti-chicken IgG recognizes chicken IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Chicken IgG whole molecule

Host: Goat

Reactivity: Chicken

Applications: ELISA, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against anti-Papain or anti-goat IgG F(c).

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.

Concentration: 2.0 mg/mL

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH_2PO_4 , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. Do not Freeze.