

Datasheet

CCKAR polyclonal antibody

Catalog Number: PAB26124

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of CCKAR.

Immunogen: A synthetic peptide corresponding to 16 amino acids at N-terminus of human CCKAR.

Host: Rabbit

Reactivity: Human

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: BLAST analysis of the peptide immunogen showed no homology with other human proteins.

Form: Liquid

Purification: Immunoaffinity chromatography

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (3-6.8 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 886

Gene Symbol: CCKAR

Gene Alias: CCK-A, CCK1-R, CCKRA

Gene Summary: This gene encodes a G-protein

coupled receptor that binds non-sulfated members of the cholecystokinin (CCK) family of peptide hormones. This receptor is a major physiologic mediator of pancreatic enzyme secretion and smooth muscle contraction of the gallbladder and stomach. In the central and peripheral nervous system this receptor regulates satiety and the release of beta-endorphin and dopamine. [provided by RefSeq]